


Premium Sutures





We're Securos Surgical,
the standard for
customer experience.

Who is Securos Surgical?

Securos is an innovator of practical solutions for veterinary surgeons. We specialise in orthopaedics, but further equip the surgical suite by offering quality daily-use items including surgical instruments and suture.

What makes us different? Why should you do business with us?

You can trust us for innovation, quality and value. And, as important as the products we bring to market, we deliver these products with the industry's most personalised service.

- Customer-focused in every aspect of what we do and every decision we make
- Problem-solving and inviting of collaboration with our customers to solve those problems
- Experts in our domain, but humble enough to constantly seek feedback
- Clear in our commitment to practical innovation, quality, value and service
- Rooted in the day-to-day reality of veterinary surgeons and general practitioners

Orthopaedics

Developed from the insights of practicing veterinary surgeons.

Instruments

Performance you can feel the moment each instrument reaches your hand.

Sutures

Premium quality at price points that drive your success.

Table of Contents

Suture Profiles	4-5
Quality and Testing	6-7
Suture Needles	8-10
Securos 300™ Needles	
Needle Types	
Cobra Black® Needles	
RiverClose™ Skin Stapler.....	11
Packaging Design and Comparison	12-13
Absorbable Suture	14-19
Securosorb™ Quick (PGA)	
Securosorb™ (PGLA)	
Securocryl™ (PGCL)	
Securdox™ (PDO)	
Antibacterial “PLUS” Absorbable Suture	20-23
Securosorb™ PLUS (PGA)	
Securdox™ PLUS (PDO)	
Securocryl™ PLUS (PGCL)	
Non-Absorbable Suture	24-25
Securolon™ (NYLON)	
Securopro™ (POLYPROPYLENE)	
Visipro (POLYPROPYLENE)	
Suture Cassettes.....	26-27
Securosorb™ (PGA)	
Securocryl™ (PGCL)	
Securodox™ (PDO)	
Chromic Gut	
Visipro™ (POLYPROPYLENE)	
Securomid™ (NYLON)	
Comparison and Cross Reference	28-31
Material Comparison	
SKU Quick Reference Charts	



Suture Profiles

Absorbable Suture



SECUROSORB™ QUICK

Fast-absorbing, Coated, Braided, Polyglycolic Acid (PGA), Synthetic Absorbable Suture

Tensile: 60% up to 7 days

Complete Absorption: 40-60 days

Indications: Skin, Oral, Ophthalmic-wedge resection

Undyed



SECURODOX™

Monofilament, Polydioxanone (PDO) Synthetic Absorbable Suture

Tensile: 75-95% up to 14 days

Complete Absorption: 180-240 days

Indications: Fascia and Ligaments, Orthopaedics-Joint Capsule, Tendon Repair, Abdominal Closure - Linea Alba, Gastric-intestinal, Vascular

Violet



SECUROSORB™

Coated, Braided, Polyglactin 910 (PGLA), Synthetic Absorbable Suture

Tensile: 65-80% up to 14 days

Complete Absorption: 60-90 days

Indications: Subcutaneous and Intradermal Suturing, Gastro-intestinal, Reproductive Tract, Muscles

Violet or Undyed



CHROMIC GUT™

Chromic Treated, Natural Absorbable Suture (Cassette only)

Tensile: Variable

Complete Absorption: Variable

Indication: Soft tissue approximation

Beige



SECUROCRYL™

Monofilament, Poly (glycolide-co-caprolactone) (PGCL) Synthetic Absorbable Suture

Tensile: 25-40% up to 14 days

Complete Absorption: 80-120 days

Indications: Subcutaneous and Intracutaneous Suturing, Gastro-intestinal, Urological and Reproductive Track, Oral Cavity Mucosa

Violet

Non-Absorbable Suture



SECUROLON™

Monofilament Nylon Non-Absorbable Suture

Indications: Skin closure, cardiovascular, ophthalmic, and neurological procedures.

Black or Blue



SECUROMID™

Braided, Cable, Nylon Non-Absorbable Suture

Indications: Skin closure

White



SECUROPRO™

Monofilament Polypropylene Non-Absorbable Suture

Indications: Skin Closure, Cardiovascular, Ophthalmology, Hernia Repairs, Tendon repair

Blue



VISIPRO™

Monofilament Polypropylene Fluorescent Pink Non-Absorbable Suture (Cassette only)

Indications: Skin Closure

Pink

Antibacterial Absorbable Suture



SECUROSORB™ PLUS

ANTIBACTERIAL Coated, Braided,
Polyglactin 910 (PGLA), Synthetic
Absorbable Suture

Tensile: 70-80% up to 14 days
Complete Absorption: 60-90 days
Indications: Subcutaneous and Intradermal
Suturing, Gastro-intestinal, Reproductive Tract,
Muscles

Violet



SECUROCRYL™ PLUS

ANTIBACTERIAL Monofilament, Poly
(glycolide-co-caprolactone) (PGCL)
Synthetic Absorbable Suture

Tensile: 25-40% up to 14 days
Complete Absorption: 80-120 days
Indications: Subcutaneous and Intracutaneous
Suturing, Gastro-intestinal, Urological and
Reproductive Track, Oral Cavity Mucosa

Violet




SECURODOX™ PLUS

ANTIBACTERIAL

Monofilament, Polydioxanone (PDO)
Synthetic Absorbable Suture

Tensile: 75-95% up to 14 days
Complete Absorption: 180-240 days
Indications: Fascia and Ligaments,
Orthopaedics-Joint Capsule, Tendon Repair,
Abdominal Closure - Linea Alba, Gastric-
intestinal, Vascular

Violet



Premium quality. Uncompromising value.



Environmentally
Controlled Clean Rooms

Regulatory Compliance

How do we build a great suture?

Quality systems at partnered manufacturer in USA.

ISO 13485 (2016): Medical Devices - Quality Management Systems

CFR Title 21 Part 820 ▪ 93/42/EEC: Medical Device Directive ▪ SOR/98 ▪ RDC 16/2013

United States Pharmacopeia ▪ 70+ Countries ▪ FDA Registered ▪ 510k(s) ▪ CE Mark (Pending)

Inspection

How do we know?

Process Verification

- 100% In-Process Needle Attachment

Sterilization Release

- Parametric
- Biological

Release Testing per United States Pharmacopeia (USP) Specifications

- (5 pcs) Needle Attachment Sutures-Needles <871>
- (10 pcs) Tensile Strength Tensile Strength <881>
- (10 pcs) Diameter/Length Testing Suture-Diameter <861>
- Visual Inspection
- Pouch Seal Verification >1.3lbf

Process Trending

- Weekly QA
- Users
- Feedback
- Clinical




Suture Performance


How do we compare?

USP<881> is the industry standard for testing tensile strength of surgical suture. The Securos brand PGCL, Securocryl, outperforms USP requirements and the market leader brand.

USP <881> Knot-Pull Tensile 2/0 PGCL Suture




USP <881> Knot-Pull Tensile 3/0 PGCL Suture




Quality Matters

That's why every single one of our suture options is made in the United States with premium materials, industry-leading tensile strength, and needles that require 12-20% less force to push through tissue than comparable products on the market. We only make premium products, but don't charge premium prices. We support your hospital's financial health, and the health of your patients, without sacrificing quality or performance.



Tensile Strength Comparison Chart

- SECUROSORB™ QUICK
- SECUROSORB™
- SECURODOX™
- SECUROCRYL™
- CHROMIC GUT



COMPLETE ABSORPTION

- 40-60 Days
- 60-90 Days
- 180-240 Days
- 80-120 Days
- Variable

Suture Needles

Securos 300™ Premium Suture Needles

Our needles are made from a propriety steel blend which is mechanically processed to make the material harder while still retaining high tensile strength. This enables the needles to have very sharp points and great durability while maintaining their flexibility. The same steel blend in our needles is used in the demanding human cardiac and plastic surgery areas where it has displaced more conventional needles.

Other leading brands reserve 300 series stainless steel needles for their higher quality lines, such as human cardiovascular. Because we believe our customers deserve premium quality, Securos offers 300 series stainless steel needles with every suture.




Needle Sharpness & Strength

Our needles are made from 300 series stainless steel, which has greater tensile strength, higher ductility and less brittleness than 400 series stainless steel.

STRENGTH

100% increase in ductility


Less likely to break when reshaping



SHARPNESS

12% less force needed

Less force, less tissue trauma




Data on File


Needle Types

REVERSE CUTTING ▽

3/8 Circle



1/2 Circle



CONVENTIONAL CUTTING ▲

1/2 Circle




SPATULA ▽

3/8 Circle




PREMIUM REVERSE CUTTING ▽

3/8 Circle




TAPERSHARP™ ▽

1/2 Circle



TAPER POINT ⊙

1/2 Circle




TAPER POINT ⊙

3/8 Circle




Cobra Black® Needles



The Cobra Black® needle is a patent pending design that enhances visibility especially in blood and smaller surgery sites. The Cobra Black® needles eliminate the hindering glare from the operating room lighting. The needles are more distinguishable than the standard stainless steel needles. Cobra Black® needles are featured on the Securosorb Quick™ suture commonly used in dental procedures, with several other selections as well.

Cobra Black® vs. Standard Needles



RiverClose™ Skin Stapler

Stainless Steel Skin Stapler

Our RiverClose™ sterile disposable skin stapler is ready-to-use for quick and easy skin closures.

RiverClose™ Features

COMFORTABLE SQUEEZE HANDLE DESIGN


READY-TO-USE STERILE DISPOSABLE STAPLER

SINGLE USE / NON-RELOADABLE

CONTAINS 35 / 316L STAINLESS STEEL STAPLES IN EITHER THE REGULAR OR WIDE STAPLE DESIGN

GOOD VISIBILITY ALLOWING FOR ACCURATE STAPLE PLACEMENT

VERY ECONOMICAL



For Veterinary Use Only

UK CODE	PRODUCT CODE	DESCRIPTION	STAPLE SIZE	PACKAGED
30494602	RP35R	35 Regular	5.4mm x 3.6mm	1 Stapler
30494613	RP35W	35 Wide	7.0mm x 4mm	1 Stapler

Packaging Design

Features and Benefits



Elegant Suture Box

- Colour coded for easy identification.
- Easy open box over-wrap pull tab offers swift and painless opening.
- Twelve pack size is optimal for storage after use.

COLOUR IDENTIFICATION

POUCH PULL SLOT


Convenient Suture Peel Pouch

- Quality Every Time. Foil pouch offers complete moisture and UV barrier.
- Easy to open Securos peel design.
- Single foil pouch unit. Reduces waste and steps to get product onto sterile field.



PROTECTIVE FOIL POUCH

EASY PEEL DESIGN



Premium Suture Racetrack

- High Performance Racetrack card significantly reduces suture memory.
- Optimal Handling. Racetrack alleviates the hassle involved with tugging and straightening the suture.
- Needle Quality. Premium needle tip protected in package by foam.

SHARPER NEEDLES

LESS MEMORY

OVAL RACETRACK CARD

Packaging Comparison

Securos' packaging is specifically designed to make it a quick and smooth conversion from the leading brand to our product.




LEADING BRAND




SECUROS

PRODUCT CODE
SUTURE SIZE
SUTURE NEEDLE
MATERIAL



PRODUCT CODE
SUTURE SIZE
SUTURE NEEDLE
MATERIAL



Absorbable Suture

SECTION




Chart Example - Needle Type

Suture Type

					Suture Size (USP/Metric)			
Needle	Material	Colour	Length		5-0 1	4-0 1.5	3-0 2	2-0 3
XX-XX XXmm X/X circle	XXXXX	XXXXX	XXcm			SKU ♦		

♦ Cobra Black® Needle

Premium Reverse Cutting ▽

SECUROSORB™ QUICK (PGA)

					Suture Size (USP/Metric)			
Needle	Material	Colour	Length		5-0 1	4-0 1.5	3-0 2	2-0 3
P-1 • 11mm 3/8 circle	Braided	Undyed	45cm		VGQ9915C ♦			

Reverse Cutting ▽

SECUROSORB™ QUICK (PGA)

					Suture Size (USP/Metric)			
Needle	Material	Colour	Length		5-0 1	4-0 1.5	3-0 2	2-0 3
FS-2 • 19mm 3/8 circle	Braided	Undyed	75cm			VGQ99330C ♦		

Conventional Cutting △

SECUROSORB™ QUICK (PGA)

					Suture Size (USP/Metric)			
Needle	Material	Colour	Length		5-0 1	4-0 1.5	3-0 2	2-0 3
FS-3 Conv • 16mm 3/8 circle	Braided	Undyed	75cm				VGQ9925C ♦	
FS-2 Conv ♦ 19mm 3/8 circle	Braided	Undyed	75cm				VGQ9931C ♦	

♦ Cobra Black® Needle

Spatula

SECUROSORB™ (PGLA)

				Suture Size (USP/Metric)	
Needle	Material	Colour	Length	8-0 0.4	7-0 0.5
Micro Point Double TG140-8 6.5mm 3/8 circle 	Braided	Violet	30cm	V9560	

Premium Reverse Cutting

SECUROSORB™ (PGLA)

				Suture Size (USP/Metric)		
Needle	Material	Colour	Length	8-0 0.4	7-0 0.5	6-0 0.7
P-1 11mm 3/8 circle 	Braided	Undyed	45cm			V9500

Conventional Cutting

SECUROSORB™ (PGLA)

				Suture Size (USP/Metric)			
Needle	Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3
FS-3 Conv 16mm 3/8 circle 	Braided	Violet	75cm	V9442			
FS-3 Conv 16mm 3/8 circle 	Braided	Violet	75cm		V9443		
FS-3 Conv 16mm 3/8 circle 	Braided	Violet	75cm			V9444	

Reverse Cutting

SECUROSORB™ (PGLA)

				Suture Size (USP/Metric)						
Needle	Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
FS-2 19mm 3/8 circle 	Braided	Violet	45cm		V9386					
FS 26mm 3/8 circle 	Braided	Violet	75cm			V9890				
FSLX 36mm 3/8 circle 	Braided	Violet	75cm				V9390			
CP 40mm 1/2 circle 	Braided	Violet	75cm						V9321	

TaperSharp™

SECUROSORB™ (PGLA)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
V-7 26mm 1/2 circle		Braided	Violet	75cm			V9180				
V-7 26mm 1/2 circle		Braided	Violet	75cm				V9181			
V-30 31mm 1/2 circle		Braided	Violet	75cm				V9360			
V-34 36mm 1/2 circle		Braided	Violet	75cm				V9363			
V-30 31mm 1/2 circle		Braided	Violet	75cm					V9361		
V-30 31mm 1/2 circle		Braided	Violet	75cm						V9362	

Taper Point

SECUROSORB™ (PGLA)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
SH-2 Plus 20mm 1/2 circle		Braided	Violet	75cm			V9114				
SH 26mm 1/2 circle		Braided	Violet	75cm				V9121			
MH-1 Plus 31mm 1/2 circle		Braided	Violet	75cm					V9138		

Non-Needle

SECUROSORB™ (PGLA)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
None		Braided	Violet	1.5m				V9025	V9026	V9027	V9028

Reverse Cutting ▽

SECUROCRYL™ (PGCL)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
P-3 13mm 3/8 circle		Monofilament	Violet	45cm	Q3203						
PS-3 16mm 3/8 circle		Monofilament	Violet	45cm		Q3205					
FS-2 19mm 3/8 circle		Monofilament	Violet	75cm	Q3209						
FS-2 19mm 3/8 circle		Monofilament	Violet	90cm		Q3201					
FS-2 19mm 3/8 circle		Monofilament	Violet	90cm			Q3202				
FS 26mm 3/8 circle		Monofilament	Violet	90cm			Q3326				
FS 26mm 3/8 circle		Monofilament	Violet	75cm				Q3327			

TaperSharp™

SECUROCRYL™ (PGCL)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
V-7 26mm 1/2 circle		Monofilament	Violet	75cm			Q3439				
V-7 26mm 1/2 circle		Monofilament	Violet	75cm				Q3440			
V-34 36mm 1/2 circle		Monofilament	Violet	75cm					Q3489		

Taper Point

SECUROCRYL™ (PGCL)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
RB-1 17mm 1/2 circle		Monofilament	Violet	75cm		Q3435					
BB 18mm 3/8 circle		Monofilament	Violet	75cm		Q3548					
RB-1 17mm 1/2 circle		Monofilament	Violet	75cm			Q3437				

Reverse Cutting ▽

SECURODOX™ (PDO)

					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS 26mm 3/8 circle		Monofilament	Violet	75cm			D9957			
PS 26mm 3/8 circle		Monofilament	Violet	75cm				D9626		

TaperSharp™ Ⓢ

SECURODOX™ (PDO)

					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
V-34 36mm 1/2 circle		Monofilament	Violet	75cm				D9380		
V-34 36mm 1/2 circle		Monofilament	Violet	90cm					D9381	
V-37 40mm 1/2 circle		Monofilament	Violet	90cm					D9366	


Taper Point Plus Ⓢ

SECURODOX™ (PDO)

					Suture Size (USP/Metric)						
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
SH-2 Plus 20mm 1/2 circle		Monofilament	Violet	75cm		D9115					
SH-2 Plus 20mm 1/2 circle		Monofilament	Violet	75cm			D9116				
SH Plus 26mm 1/2 circle		Monofilament	Violet	75cm			D9124				
SH Plus 26mm 1/2 circle		Monofilament	Violet	75cm				D9125			
MH-1 Plus 31mm 1/2 circle		Monofilament	Violet	75cm				D9133			

Antibacterial Absorbable PLUS Suture

SECTION




Antibacterial Suture

Plus Surgical Suture Products

- ✚ Securos PLUS sutures provide the same handling and performance of our other suture products, but with the added benefit of antibacterial protection.
- ✚ Securos PLUS sutures are coated with chlorhexidine diacetate, a proven antibacterial agent, that result in a zone of inhibition around the suture line.
- ✚ Our PLUS antibacterial agent is recognised as being effective against bacteria such as Staphylococcus aureus, Staphylococcus epidermidis, MRSA and MRSE.

A great alternative
for high-risk patients
or procedures where
bacteria is a concern.

Veterinary Use Only



SECUROSORB™ PLUS
ANTIBACTERIAL Coated,
Braided, Polyglactin 910 (PGLA),
Synthetic Absorbable Suture

Violet



SECURODOX™ PLUS
ANTIBACTERIAL Monofilament,
Polydioxanone (PDO) Synthetic
Absorbable Suture

Violet



SECUROCRYL™ PLUS
ANTIBACTERIAL Monofilament,
Poly (glycolide-co-caprolactone)
(PGCL) Synthetic Absorbable
Suture

Violet

Zone of Inhibition




PLUS Suture product with
Antibacterial coating

* Petri dish images is for illustrative
purposes only, zone of inhibition
testing results can vary



Suture product without
Antibacterial coating

* Petri dish images is for illustrative
purposes only, zone of inhibition
testing results can vary




MRSE - Violet Plus Suture: size 2/0



S. aureus- Violet Plus Suture: size 2/0

TENSILE STRENGTH CHART



Images of Securos PLUS Suture product creating a zone of inhibition in vitro, reducing Staph aureus and MRSE bacterial colonization.

Reverse Cutting ▽

SECUROSORB™ Plus

Antibacterial					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-2 19mm 3/8 circle		Braided	Violet	75cm		VP422				
FS-1 24mm 3/8 circle		Braided	Undyed	75cm			VP442			
CP-1 37mm 3/8 circle		Braided	Violet	75cm				VP466		
CP-1 37mm 3/8 circle		Braided	Undyed	75cm					VP267	

Taper Point ⊙

SECUROSORB™ Plus

Antibacterial					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
SH 26mm 1/2 circle		Braided	Violet	75cm			VP316			

Premium Reverse Cutting ▽

SECUROCRYL™ Plus

Antibacterial					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-2 19mm 3/8 circle		Braided	Violet	75cm		QP922				

Taper Point ⊙

SECUROCRYL™ Plus

Antibacterial					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
SH 26mm 1/2 circle		Braided	Violet	75cm			QP316			
SH 26mm 1/2 circle		Braided	Violet	75cm				QP317		
CT-1 36mm 1/2 circle		Braided	Violet	75cm					QP346	

Reverse Cutting ▽

SECURODOX™ Plus

Antibacterial					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-2 24mm 3/8 circle		Braided	Violet	75cm			DP452			
CP-2 26mm 3/8 circle		Braided	Violet	75cm				DP969		
CP-1 37mm 3/8 circle		Braided	Violet	75cm					DP467	

Non- Absorbable Suture

SECTION



Reverse Cutting ▽

SECUROLON™ (NYLON)

					Suture Size (USP/Metric)					
Needle		Material	Colour	Length	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-2 19mm 3/8 circle		Monofilament	Blue	75cm		662BL				
PS 26mm Prime 3/8 circle		Monofilament	Blue	75cm			1685BL			
FS 26mm 3/8 circle		Monofilament	Blue	45cm			320BL			
PS 26mm Prime 3/8 circle		Monofilament	Blue	75cm				1626BL		

Conventional Cutting △

SECUROPRO™ (POLYPROPYLENE)

					Suture Size (USP/Metric)				
Needle		Material	Colour	Length	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS-3 CONV 16mm 3/8 circle		Monofilament	Blue	75cm		P8549			

Reverse Cutting ▽

SECUROPRO™ (POLYPROPYLENE)

					Suture Size (USP/Metric)				
Needle		Material	Colour	Length	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
FS 26mm Prime 3/8 circle		Monofilament	Blue	45cm		P8684			

TaperSharp™

SECUROPRO™ (POLYPROPYLENE)

					Suture Size (USP/Metric)				
Needle		Material	Colour	Length	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
V7 DBL 26mm 1/2 circle		Monofilament	Blue	90cm			P8977		

Suture Cassettes

SECTION



Suture Micro-Cassettes

Cassette Features

EASY-OPEN LID

CONVENIENT BUILT-IN CUTTER

ERGONOMICAL DESIGN

POCKET SIZE

ECO-FRIENDLY

(50% less plastic than traditional cassette packaging)



Cassette Rack



SKU	DESCRIPTION
MCR-100	Cassette Rack

Fits 5 Suture Cassettes

SECUROSORB™	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCV-4/0	Braided	Violet	4-0 (1.5)	25m
MCV-3/0	Braided	Violet	3-0 (2)	25m
MCV-2/0	Braided	Violet	2-0 (3)	25m

Once open dispose of after 60 days

CHROMIC GUT	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCG-3/0	Twisted	Beige	3-0 (3)	75m
MCG-2/0	Twisted	Beige	2-0 (3.5)	50m
MCG-1	Twisted	Beige	1 (5)	50m
MCG-2	Twisted	Beige	2 (6)	45m
MCG-3	Twisted	Beige	3 (7)	45m

SECURODOX™	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCD-3/0	Monofilament	Violet	3-0 (2)	25m
MCD-2/0	Monofilament	Violet	2-0 (3)	25m
MCD-1	Monofilament	Violet	1 (4)	25m

Once open dispose of after 60 days

VISIPRO™	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCF-3/0	Monofilament	Pink	3-0 (2)	110m
MCF-2/0	Monofilament	Pink	2-0 (3)	110m

SECUROCRYL™	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCQ-3/0	Monofilament	Violet	3-0 (2)	25m
MCQ-2/0	Monofilament	Violet	2-0 (3)	25m
MCQ-0	Monofilament	Violet	0 (3.5)	25m

Once open dispose of after 30 days

SECUROMID™	SUTURE MATERIAL	SUTURE COLOUR	SUTURE SIZE USP (metric)	SUTURE LENGTH
MCS-3/0	Braided	White	3-0 (2)	100m
MCS-2/0	Braided	White	2-0 (3)	100m
MCS-1	Braided	White	1 (4)	75m
MCS-3	Braided	White	3 (6)	50m

Comparison and Cross- Reference

SECTION



Material Comparison Charts

Absorbable Suture

SECUROS	ETHICON	COVIDIEN
SECUROSORB™ QUICK	Vicryl Rapide™*	Velosorb™* Fast Suture
SECUROSORB™	Vicryl™*	Polysorb™* Suture
SECUROCRYL™	Monocryl™*	Biosyn™* Suture
SECURODOX™	PDS™* II	Maxon™* Suture
CHROMIC GUT	Chromic Catgut	Chromic Catgut
SECUROSORB™ PLUS	Vicryl™* Plus	—
SECUROCRYL™ PLUS	Monocryl™* Plus	—
SECURODOX™ PLUS	PDS™* Plus	—

Non-Absorbable Suture

SECUROS	ETHICON	COVIDIEN
SECUROLON™	Ethilon™*	Monosof™* Suture Dermalon™* Suture
SECUROPRO™	Prolene™*	Surgipro™* Suture Surgipro™* II Suture
VISIPRO™	—	—
SECUROMID™	Nurolon™*	Surgilon™* Suture

™* Third party brands are trademarks of their respective owners.

Quick Reference Charts—Ordering

Absorbable Suture

SECUROSORB™ (PGLA)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494271	V9560	0.4 (8/0)	Micro Point Double 3/8 Circle 6.5mm	30cm
30494282	V9500	0.7 (6/0)	Premium Reverse Cutting 3/8 Circle 11mm	45cm
30494293	V9442	1 (5/0)	Conventional Cutting 3/8 Circle 16mm	75cm
30494307	V9443	1.5 (4/0)	Conventional Cutting 3/8 Circle 16mm	75cm
30494318	V9386	1.5 (4/0)	Reverse Cutting 3/8 Circle 19mm	45cm
30494329	V9114	2 (3/0)	Taper Point + 1/2 Circle 20mm	75cm
30494330	V9180	2 (3/0)	Tapersharp 1/2 Circle 26mm	75cm
30494340	V9444	2 (3/0)	Conventional Cutting 3/8 Circle 16mm	75cm
30494351	V9890	2 (3/0)	Reverse Cutting 3/8 Circle 26mm	75cm
30494362	V9121	3 (2/0)	Taper Point + 1/2 Circle 26mm	75cm
30494384	V9025	3 (2/0)	Non Needled	1.5m
30494395	V9181	3 (2/0)	Tapersharp 1/2 Circle 26mm	75cm
30494409	V9360	3 (2/0)	Tapersharp 1/2 Circle 31mm	75cm
30494410	V9363	3 (2/0)	Tapersharp 1/2 Circle 36mm	75cm
30494420	V9390	3 (2/0)	Reverse Cutting 3/8 Circle 36mm	75cm
30494431	V9361	3.5 (0)	Tapersharp 1/2 Circle 31mm	75cm
30494373	V9138	3.5 (0)	Taper Point + 1/2 Circle 31mm	75cm
30494442	V9026	3.5 (0)	Non Needled	1.5m
30494453	V9321	4 (1)	Reverse Cutting 1/2 Circle 40mm	75cm
30494464	V9362	4 (1)	Tapersharp 1/2 Circle 31mm	75cm
30494475	V9027	4 (1)	Non Needled	1.5m
30494486	V9028	5 (2)	Non Needled	1.5m
SECUROSORB™ Quick (PGA)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494497	VGQ9915C	1 5/0	Premium Reverse Cutting 3/8 Circle 11mm	45cm
30494500	VGQ9930C	1.5 4/0	Reverse Cutting 3/8 Circle 19mm	75cm
30494511	VGQ9925C	2 3/0	Conventional Cutting 3/8 Circle 16mm	75cm
30494522	VGQ9931C	2 (3/0)	Conventional Cutting 3/8 Circle 19mm	75cm
SECUROCRYL™ (PGCL)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494045	Q3203	1 (5/0)	Reverse Cutting 3/8 Circle 13mm	45cm
30494056	Q3209	1 (5/0)	Reverse Cutting 3/8 Circle 19mm	75cm
30494067	Q3435	1.5 (4/0)	Taper Point 1/2 Circle 17mm	75cm
30494078	Q3548	1.5 (4/0)	Taper Point 3/8 Circle 17mm	75cm
30494089	Q3205	1.5 (4/0)	Reverse Cutting 3/8 Circle 16mm	45cm
30494090	Q3201	1.5 4/0	Reverse Cutting 3/8 Circle 19mm	90cm
30494103	Q3202	2 (3/0)	Reverse Cutting 3/8 Circle 19mm	90cm
30494114	Q3326	2 (3/0)	Reverse Cutting 3/8 Circle 26mm	90cm
30494125	Q3437	2 (3/0)	Taper Point 1/2 Circle 17mm	75cm
30494136	Q3439	2 (3/0)	Tapersharp 1/2 Circle 26mm	75cm
30494147	Q3327	3 (2/0)	Reverse Cutting 3/8 Circle 26mm	75cm
30494158	Q3440	3 (2/0)	Tapersharp 1/2 Circle 26mm	75cm
30494169	Q3489	3.5 (0)	Tapersharp 1/2 Circle 36mm	75cm

SECURODOX™ (PDO)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494170	D9115	1.5 (4/0)	Taper Point + 1/2 Circle 20mm	75cm
30494180	D9116	2 (3/0)	Taper Point + 1/2 Circle 20mm	75cm
30494191	D9124	2 (3/0)	Taper Point + 1/2 Circle 26mm	75cm
30494205	D9957	2 (3/0)	Reverse Cutting 3/8 Circle 26mm	75cm
30494216	D9125	3 (2/0)	Taper Point + 1/2 Circle 26mm	75cm
30494227	D9133	3 (2/0)	Taper Point + 1/2 Circle 31mm	75cm
30494238	D9626	3 (2/0)	Reverse Cutting 3/8 Circle 26mm	75cm
30494249	D9380	3 (2/0)	Tapersharp 1/2 Circle 36mm	75cm
30494250	D9381	3.5 (0)	Tapersharp 1/2 Circle 36mm	90cm
30494260	D9366	3.5 (0)	Tapersharp 1/2 Circle 40cm	90cm
SECUROSORB™ PLUS				
ANTIBACTERIAL				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30476940	VP422	1.5 (4/0)	Reverse Cutting 3/8 Circle 19mm	75cm
30476950	VP442	2 (3/0)	Reverse Cutting 3/8 Circle 24mm	75cm
30476881	VP316	2 (3/0)	Taper Point 1/2 Circle 26mm	75cm
30476961	VP466	3 (2/0)	Reverse Cutting 3/8 Circle 36mm	75cm
30476972	VP267	3.5 (0)	Reverse Cutting 3/8 Circle 36mm	75cm
SECUROCRYL™ Plus				
ANTIBACTERIAL				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30476870	QP922	1.5 (4/0)	Reverse Cutting 3/8 Circle 19mm	75cm
30476881	QP316	2 (3/0)	Taper Point 1/2 Circle 26mm	75cm
30476892	QP317	3 (2/0)	Taper Point 1/2 Circle 26mm	75cm
30476906	QP346	3.5 (0)	Taper Point 1/2 Circle 36mm	75cm
SECURODOX™ Plus				
ANTIBACTERIAL				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30476917	DP452	2 (3/0)	Reverse Cutting 3/8 Circle 24mm	75cm
30476928	DP969	3 (2/0)	Reverse Cutting 3/8 Circle 26mm	75cm
30476939	DP467	3.5 (0)	Reverse Cutting 3/8 Circle 36mm	75cm

Non-Absorbable Suture

SECUROLON™ (NYLON)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494533	662BL	1.5 (4/0)	Reverse Cutting 3/8 Circle 19mm	75cm
30494544	1685BL	2 (3/0)	Reverse Cutting Prime 3/8 Circle 26mm	75cm
30494555	320BL	2 (3/0)	Reverse Cutting 3/8 Circle 26mm	45cm
30494566	1626BL	3 (2/0)	Reverse Cutting Prime 3/8 Circle 26mm	75cm
SECUOPRO™ (POLYPROPYLENE)				
MWI	SECUROS	SIZE	NEEDLE	LENGTH
30494588	P8549	2 (3/0)	Conventional Cutting 3/8 Circle 16mm	75cm
30494599	P8684	2 (3/0)	Reverse Cutting 3/8 Circle Prime 26mm	45cm
30494577	P8977	3 (2/0)	Tapersharp 1/2 Circle 26mm Double	90cm

Quick Reference Charts—Comparison

Absorbable Suture

MWI	Securosorb™	Ethicon Vicryl™
30494271	V9560	W9560
30494282	V9500	W9500T
30494293	V9442	W9442
30494307	V9443	W9443
30494318	V9386	W9386
30494329	V9114	W9114
30494330	V9180	W9180
30494340	V9444	W9444
30494351	V9890	W9890
30494362	V9121	W9121
30494384	V9025	W9025
30494395	V9181	W9181
30494409	V9360	W9360
30494410	V9363	W9363
30494420	V9390	W9390
30494431	V9361	W9361
30494373	V9138	W9138
30494442	V9026	W9026
30494453	V9321	W9321
30494464	V9362	W9362
30494475	V9027	W9027
30494486	V9028	W9028
30494497	VGQ9915C	W9915
30494500	VGQ9930C	W9930
30494511	VGQ9925C	W9925
30494522	VGQ9931C	W9931
MWI	Securocryl™	Ethicon Monocryl™
30494045	Q3203	W3203
30494056	Q3209	W3209
30494067	Q3435	W3435
30494078	Q3548	W3548
30494089	Q3205	W3205
30494090	Q3201	W3201
30494103	Q3202	W3202
30494114	Q3326	W3326
30494125	Q3437	W3437
30494136	Q3439	W3439
30494147	Q3327	W3327
30494158	Q3440	W3440
30494169	Q3489	W3489
MWI	Securodox™	Ethicon PDS™ II
30494170	D9115	W9115H
30494180	D9116	W9116H
30494191	D9124	W9124H
30494205	D9957	W9957T
30494216	D9125	W9125H
30494227	D9133	W9133H
30494238	D9626	W9626T
30494249	D9380	W9380H
30494250	D9381	W9381H
30494260	D9366	W9366T

Non-Absorbable Suture

MWI	Securodox™	Ethicon Ethilon™
30494533	662BL	W319
30494544	1685BL	W1685T
30494555	320BL	W320
30494566	1626BL	W1626T
MWI	Securodox™	Ethicon Prolene™
30494577	P8977	W8977
30494588	P8549	W8549
30494599	P8684	W8684
MWI	Securodox™	Ethicon Proximate™
30494602	RP35R	PMW35
30494613	RP35W	PMR35



Where knowledge,
reach and partnership
shape healthcare delivery.